

Dulbecco's Modified Eagle Medium, Low Glucose

DMEM 低糖液体培养基（含双抗,含丙酮酸钠,不含 HEPES）

产品编号	产品名称	规格
BL1916A	DMEM低糖液体培养基（含双抗,含丙酮酸钠,不含HEPES）	500ml

产品简介:

本产品含青霉素和链霉素，并经 0.22μm 滤膜过滤除菌，广泛应用于细胞的培养。

产品组成:

Components	Molecular Weight	Concentration (mg/L)	mM
Amino Acids			
Glycine	75.0	30.0	0.4
L-Arginine hydrochloride	211.0	84.0	0.39810428
L-Cystine 2HCl	313.0	63.0	0.20127796
L-Glutamine	146.0	584.0	4.0
L-Histidine hydrochloride-H ₂ O	210.0	42.0	0.2
L-Isoleucine	131.0	105.0	0.8015267
L-Leucine	131.0	105.0	0.8015267
L-Lysine hydrochloride	183.0	146.0	0.7978142
L-Methionine	149.0	30.0	0.20134228
L-Phenylalanine	165.0	66.0	0.4
L-Serine	105.0	42.0	0.4
L-Threonine	119.0	95.0	0.79831934
L-Tryptophan	204.0	16.0	0.078431375
L-Tyrosine disodium salt dihydrate	261.0	104.0	0.39846742
L-Valine	117.0	94.0	0.8034188
Vitamins			
Choline chloride	140.0	4.0	0.028571429
D-Calcium pantothenate	477.0	4.0	0.008385744
Folic Acid	441.0	4.0	0.009070295
Niacinamide	122.0	4.0	0.032786883
Pyridoxine hydrochloride	204.0	4.0	0.019607844
Riboflavin	376.0	0.4	0.0010638298
Thiamine hydrochloride	337.0	4.0	0.011869436
i-Inositol	180.0	7.2	0.04

Note: For in vitro research use only, not for diagnostic or therapeutic use, This product is not a medical device.
 注意: 在体外研究使用, 不用于诊断或治疗用途, 本产品不是医疗装置。



Inorganic Salts

Calcium Chloride (CaCl ₂) (anhyd.)	111.0	200.0	1.8018018
Ferric Nitrate (Fe(NO ₃) ₃ ·9H ₂ O)	404.0	0.1	2.4752476E-4
Magnesium Sulfate (MgSO ₄) (anhyd.)	120.0	97.67	0.8139166
Potassium Chloride (KCl)	75.0	400.0	5.3333335
Sodium Bicarbonate (NaHCO ₃)	84.0	3700.0	44.04762
Sodium Chloride (NaCl)	58.0	6400.0	110.344826
Sodium Phosphate monobasic (NaH ₂ PO ₄ ·H ₂ O)	138.0	125.0	0.9057971

Other Components

D-Glucose (Dextrose)	180.0	1000.0	5.5555553
Phenol Red	376.4	15.0	0.039851222
Sodium Pyruvate	110.0	110.0	1.0

注意事项:

- 1、本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 2、为了您的安全和健康，请穿实验服并戴一次性手套操作。

保存条件:

2-8℃，避光保存，有效期 12 个月。

